

## Bind-It™

### Designed for Nuclear Medicine

#### The Next Generation of Radioactive Decontaminants

Reasons to use Bind-It Decontaminant Hand Soap every time you wash your hands:

- Will not dry your hands like ordinary lab soap
- Fresh, clean scent
- Economical & convenient
- Touch-free wall mount dispenser
- Meets the intent of ALARA
- Removes contamination every time you wash

When you use Bind-It Decontaminant Spray the difference is clear:

- Will NOT cloud lead glass
- Safely decontaminate Syringe Shields & L-Blocks
- Removes dirt, dust, oils and grease as well as Radioactive Contamination
- Used regularly it provides a protective barrier that aids in subsequent cleanup
- Very mild, pleasant scent



LTI047

LTI046

LTI030

LTI051

Bind-It Decontamination Introductory Kit

(1) Touch-free dispenser

(1) Bind-It Hand Soap- 32 oz bottle

(1) Bind-It Decontaminant Spray - 32 oz bottle

LTI046

Bind-It Decontaminant Hand Soap - 8 oz bottles, Case of 12

LTI047

Bind-It Decontaminant Hand Soap - 32 oz bottles, Case of 12

6-11

1-5



LTI030

Bind-It Decontaminant Spray Cleaner - 32 oz bottles, Case of 12

6-11

1-5

LTI031

Decontaminant Wipes, 50/Package



LTIHCP

Patient Care Pack

Includes: 1 Bottle of 16 oz S4FE-D CRBN Decon Soap  
1 Bottle of 32 oz S4FE-D CRBN Decon Spray  
1 Bottle of 16 oz S4FE-D CRBN Decon Concentrate  
1 Bag of S4FE-D CRBN Decon Wipes