

## Bind-It™

### Designed for Nuclear Medicine

#### The Next Generation of Radioactive Decontaminants

Reasons to use Bind-It Decontaminant Hand Soap every time you wash your hands:

- Will not dry your hands like ordinary lab soap
- Fresh, clean scent
- Economical & convenient
- Touch-free wall mount dispenser
- Meets the intent of ALARA
- Removes contamination every time you wash

When you use Bind-It Decontaminant Spray the difference is clear:

- Will NOT cloud lead glass
- Safely decontaminate Syringe Shields & L-Blocks
- Removes dirt, dust, oils and grease as well as Radioactive Contamination
- Used regularly it provides a protective barrier that aids in subsequent cleanup
- Very mild, pleasant scent



LT1047

LT1016

LT1030

LT1051

LT1051

Bind-It Decontamination Introductory Kit

- (1) Touch-free dispenser
- (1) Bind-It Hand Soap- 32 oz bottle
- (1) Bind-It Decontaminant Spray - 32 oz bottle

LT1016

Decontaminant Hand Soap - 16 oz bottles, Case of 12

LT1047

S4FE-D Decontaminant Hand Soap - 32 oz Pouch, Case of 12

6-11

1-5



LT1030

Bind-It Decontaminant Spray Cleaner - 32 oz bottles, Case of 12

6-11

1-5

LT1031

Decontaminant Wipes, 50/Package

Patient Care Packs allow patients to safely and effectively minimize and remove contamination from their isolation areas.

- The 16oz Hand Soap will accommodate many hand washings over the course of the patient's post treatment period and can be used in the shower
- The Ready-To-Use Spray is ready to use on the bathroom sink, bathtub/shower and floor
- The Concentrate is ideal for cleaning the sink, bathtub/shower and toilet bowl (the largest typical point of contamination)



LT1HCP

Patient Care Pack

- Includes:
- 1 Bottle of 16 oz S4FE-D CRBN Decon Soap
  - 1 Bottle of 32 oz S4FE-D CRBN Decon Spray
  - 1 Bottle of 16 oz S4FE-D CRBN Decon Concentrate
  - 1 Bag of S4FE-D CRBN Decon Wipes