

Designed for Nuclear Medicine

The Next Generation of Radioactive Decontaminants



LTI047

LTI016

LTI030

Reasons to use S4FE-D Decontaminant Hand Soap every time you wash your hands:

- Will not dry your hands like ordinary lab soap
- Fresh, clean scent
- Economical & convenient
- Touch-free wall mount dispenser
- Meets the intent of ALARA
- Removes contamination every time you wash

When you use S4FE-D Decontaminant Spray the difference is clear:

- Will NOT cloud lead glass
- Safely decontaminate Syringe Shields & L-Blocks
- Removes dirt, dust, oils and grease as well as Radioactive Contamination
- Used regularly it provides a protective barrier that aids in subsequent cleanup
- Very mild, pleasant scent



LTI051

S4FE-D Decontamination Introductory Kit

(1) Touch-free dispenser

(1) S4FE-D Hand Soap- 32 oz bottle

(1) S4FE-D Decontaminant Spray - 32 oz bottle

LTI016

Decontaminant Hand Soap - 16 oz bottles, Case of 12

LTI047

S4FE-D Decontaminant Hand Soap - 32 oz Pouch, Case of 12

6-11

1-5



LTI030

S4FE-D Decontaminant Spray Cleaner - 32 oz bottles, Case of 12

6-11

1-5

LTI031

S4FE-D Decontaminant Wipes, 50/Package

Patient Care Packs allow patients to safely and effectively minimize and remove contamination from their isolation areas.



LTIHCP

Patient Care Pack

- Includes:
- 1 Bottle of 16 oz S4FE-D CRBN Decon Soap
 - 1 Bottle of 32 oz S4FE-D CRBN Decon Spray
 - 1 Bottle of 16 oz S4FE-D CRBN Decon Concentrate
 - 1 Bag of S4FE-D CRBN Decon Wipes